


The User Instruction of Dental Anesthesia Delivery Equipment Easy-III



Manufacturer : Ziyang Smiles Technology Co.,Ltd

**Domicile: Block B2/C1, 4th Floor, Building C, No.3 xiandai Road, Ziyang City,
Sichuan China**

**Manufacturer address: Block B2/C1, 4th Floor, Building C, No.3 xiandai Road,
Ziyang City, Sichuan China**

Tel/Fax: 0086 28-26241001

Mob: 0086 13804574465

Web : www.sunshine-smiles.com

Email: sunnysu@sunshine-smiles.com

European Representative:

Name:Medwheat Tech Service GmbH.

Add:Max-Plank-Straße 4 85609 Aschheim b.München Germany

Tel/fax:0049-89-716718760/0049-89-716718761

Catalogue

I、 Name Of Product.....	3
II、 Medical Use.....	3
III、 Model、 Specification.....	3
IV、 Configuration and Composition ...	3
V、 Performance.....	5
VI、 Technical Parameters.....	6
VII、 Operation Instructions.....	8
VIII、 Maintenance and Care... ..	11
IX、 Contraindication.....	11
X、 Attention.....	11
XI、 Warnings and Tips.....	14
XII、 Date of Produce.....	15
XII、 The period of Use	15

XII、 Date of Produce: Show on the label

XIII、 The period of Use: 5 years

5. During transportation and using, the product shall be protected from violent vibration, drop, impact or extrusion.

6. If the product reaches the scrapping period, the machine and accessories shall be disposed in accordance with relevant national regulations.

7. The adapter should meet the requirements of IEC 60601-1.

8. The product can not use in Oxygen-enriched state, is not available for flammable anesthetic.

XI. Warnings and Tips:

1. Make sure the O-ring installed and complete before aspiration.

2. Lithium-ion batteries are recommended to be charged once every 3 months according to standard charging and discharging methods.

3. Each main machine must be charged with the matched power adapter. The main machine will not work when charging. Other cables and accessories may negatively affect EMC performance.

4. If the equipment is abnormal, please contact our company for repair, do not disassemble it without authorization, and provide the circuit diagram, component list, drawing note and correction rules as required by our company.

5. For your safety, do not modify the ME EQUIPMENT without authorization.

Document NO.SM/IFU-02 Rav.02

Date of draft: 01th, Dec, 2018

Date of review: 17th, Feb, 2020

Revision date: 17th, Feb, 2020

I. Name of Product:

Dental Anesthesia Delivery Equipment

II. Medical Use

This device is intended for use only in local anesthetic agents for dental applications.

It's appropriate for local anesthesia during the treatments of oral diseases, such as intraoral cure, repair, reform, small size of surgical operation and so on. The device is use with 1.8ml standard Card local anesthetics (Primacaine, Scandonest) and Disposable Dental Needle

III. Model、 Specification:

Easy-III

IV. Configuration and Composition:

1. One main machine


2. Four boosting pipes (NO.1 British system, NO.2 Metric system)

3. One base

4. One power adapter (With One USB Cable)


5. O-ring ten pieces

Main Machine



USB Sock: connect the USB power adapter(DC5V,1A)to charge or connected to the computer through the USB cable to change the music.

LED Display :



Conducted RFIEC 61000-4-6	3Vrms 0.15MHz~80MHz 6Vrms in ISM bands between 0,15MHz and 80MHz 80% AM at 1kHz (only charging mode)																									
Radiated RFIEC 61000-4-3	For EM fields: 3V/m 80MHz~2.7GHz For Proximity fields from RF wireless communications equipment:																									
	<table border="1"> <thead> <tr> <th>Freq MHz</th> <th>Test Level P: max power, d:distance, E:Immunity Level</th> </tr> </thead> <tbody> <tr> <td>385</td> <td>P=1.8W d=0.3m E=27V/m for TETRA400</td> </tr> <tr> <td>450</td> <td>P=2W d=0.3m E=28V/m for GMRS460; FRS460</td> </tr> <tr> <td>710</td> <td rowspan="3">P=0.2W d=0.3m E=9V/m for LTE Band 13, 17</td> </tr> <tr> <td>745</td> </tr> <tr> <td>780</td> </tr> <tr> <td>810</td> <td rowspan="2">P=2W d=0.3m E=28V/m for GSM800/900; TETRA800; iDEN820; CDMA850;LTE Band 5</td> </tr> <tr> <td>870</td> </tr> <tr> <td>930</td> <td rowspan="2">P=2W d=0.3m E=28V/m for GSM1800, CDMA1900; GSM1900; DECT; LTE Band 1,3,4,35;UMTS</td> </tr> <tr> <td>1720</td> </tr> <tr> <td>1845</td> <td rowspan="2">P=2W d=0.3m E=28V/m for Bluetooth, WLAN 802.11 b/g/n, RFID 2450, LTE Band 7</td> </tr> <tr> <td>1970</td> </tr> <tr> <td>2450</td> <td rowspan="3">P=0.2W d=0.3m E=9V/m for WLAN 802.11 a/n</td> </tr> <tr> <td>5240</td> </tr> <tr> <td>5500</td> </tr> <tr> <td>5785</td> <td></td> </tr> </tbody> </table>	Freq MHz	Test Level P: max power, d:distance, E:Immunity Level	385	P=1.8W d=0.3m E=27V/m for TETRA400	450	P=2W d=0.3m E=28V/m for GMRS460; FRS460	710	P=0.2W d=0.3m E=9V/m for LTE Band 13, 17	745	780	810	P=2W d=0.3m E=28V/m for GSM800/900; TETRA800; iDEN820; CDMA850;LTE Band 5	870	930	P=2W d=0.3m E=28V/m for GSM1800, CDMA1900; GSM1900; DECT; LTE Band 1,3,4,35;UMTS	1720	1845	P=2W d=0.3m E=28V/m for Bluetooth, WLAN 802.11 b/g/n, RFID 2450, LTE Band 7	1970	2450	P=0.2W d=0.3m E=9V/m for WLAN 802.11 a/n	5240	5500	5785	
	Freq MHz	Test Level P: max power, d:distance, E:Immunity Level																								
	385	P=1.8W d=0.3m E=27V/m for TETRA400																								
	450	P=2W d=0.3m E=28V/m for GMRS460; FRS460																								
	710	P=0.2W d=0.3m E=9V/m for LTE Band 13, 17																								
	745																									
	780																									
	810	P=2W d=0.3m E=28V/m for GSM800/900; TETRA800; iDEN820; CDMA850;LTE Band 5																								
	870																									
	930	P=2W d=0.3m E=28V/m for GSM1800, CDMA1900; GSM1900; DECT; LTE Band 1,3,4,35;UMTS																								
	1720																									
	1845	P=2W d=0.3m E=28V/m for Bluetooth, WLAN 802.11 b/g/n, RFID 2450, LTE Band 7																								
	1970																									
2450	P=0.2W d=0.3m E=9V/m for WLAN 802.11 a/n																									
5240																										
5500																										
5785																										

3.The storage products warehouse environment temperature for - 10 °C to + 45 °C of indoor

cleaning ,to avoid contact with corrosive substances, stay away from fire and heat source.





4.Please keep away from strong magnetic field and alternating current equipment to avoid electromagnetic interference.

Electromagnetic compatibility statement as follows: The products for Professional health care facility environment.

Standard	Test methods and levels
Radiated Emissions CISPR 11	Group 1 Class B
Conducted Emissions CISPR 11	Group 1 Class B (only charging mode)
Harmonic emissions IEC 61000-3-2	Class A, (only charging mode)
Voltage fluctuations/ flicker emissions IEC 61000-3-3	Complies, (only charging mode)
Electrostatic discharge (ESD) IEC 61000-4-2	8 kVcontact 2 kV, 4 kV, 8 kV, 15 kVair
Electrical fast transient/burst IEC 61000-4-4	□2 kV 100KHz repetition frequency (only charging mode)
Surge input AC power port IEC 61000-4-5	□1kV Line-to-Line (only charging mode)
Voltage dips IEC 61000-4-11	0% UT for 0.5 cycle At 0°45°90°135°180°225°270°and 315° 0% UT for 1 cycle and 70% UT for 25 cycles Single phase at 0° (only charging mode)
Voltage interruptions IEC 61000-4-11	0 % UT; 250/300 cycle (only charging mode)
Power frequency(50/60Hz)Magnetic fieldIEC-61000-4-8	30 A/m 50Hz / 60Hz

Button Instruction

Table 1

Picture	Button Instruction
	Inductive switch: When the inductive switch is unlocked, use the finger to block the inductive switch to start inject medicine ,and the finger leaves the inductive switch to stop pushing the medicine.
	The aspiration/inductive switch unlock button: press the button to unlock (or lock) the inductive switch, the blue work indicator lights up and enters the working state; press and hold the button for 2-5 seconds to do aspiration
	Speed/Music: Press the button one second to set the High speed/Medium speed/Low speed; Pres and hold the button for 2-5 seconds to set Mute mode、 Warning tone (beep every 0.1ml)、 Music
	Press and hold the button for 2-5 seconds to turn on /off the power, Double click the button to reset the push rod.

V、 Performance:

1.Basic performance of Dental Anesthesia Delivery Equipment:

booster tube: There is no damage when falling at a height of 1.5 meters easy to insert and pull;

push time: push time error does not exceed ±10%;

push dose: push dose 1.8ml, The error does not exceed ±10%

2.The connection situation: the booster tube and the host connection should be reliable; the button should be flexible and reliable.

3.Appearance: The outer surface of the machine should be smooth, and there should be no defects such as corrosion spots, dirt, scratches, sharp edges and burrs.

4.Safety requirements: electrical safety performance fully implements the requirements of IEC 60601-1; electromagnetic compatibility should comply with the provisions of IEC 60601-1-2.

5.Electromagnetic compatibility: The electromagnetic compatibility of the dental anesthesia delivery equipment is in accordance with the provisions of IEC 60601-1-2.

6.Environmental test: The environmental test of the dental anesthesia delivery equipment is in accordance with the requirements of IEC 60601-1 and meets the standard requirements.

7. The unique features of Easy-III

7.1 The display should show clear, no black screen, defects, display confusion and other defects.

7.2 When the battery is low, the warning light flashes with an audible alert function.

7.3 When you making injection, you can select a music playback function to play a song, or you can choose to mute or tone.

VI、 Technical Parameters:

- Security feature::

Type classification of electric shock protection: For cure: Internal electric source. For charge:


Class II

Degree classification of electric shock protection: BF type of application part

VIII、 Maintenance and Care:

1. Check the O-ring every day. If there is a worn O-ring, it should be replaced immediately. If the aspiration function is not used for a long time, it is recommended to remove the O-ring with a dental probe.

2. In order to extend the battery life, it is recommended to charge the battery when there is still one grid remaining.

3. Restore the original factory settings: If there is a continuous alarm or other fault status during use, connect the host to the power adapter or remove the battery and push it in again. After the prompt tone is heard, double click  to reset the push rod, the original factory settings will be restored.

IX、 Contraindication: Intraligamentary injection is prohibited in patients with active periodontal disease.



X、 Attention:

1.This product is used at an altitude under 3000meters above with temperature of 5-40°C,relative humidity of 30-80%.This product is only applied for the local anesthesia in the Stomatology Department. The selection of anesthetic and needle shall be purchased according to doctors' clinical requirements. Any results caused by misuse of anesthetic and needle have nothing to do with the product.

2.The machine and charger should be disinfected after use, It is recommended to spray the disinfectant onto a clean paper towel to wipe the machine. Do not spray the disinfectant directly onto the machine.

(It is recommended to copy the internal four files before replacing, and then format the storage space, then move the four files into the storage space one by one)

NOTE: 001key.mp3 and 002prompt.mp3 MUST BE RETAINED !

9. **Aspiration** : Press the  to unlock the button and the blue work indicator lights up to enter working state when you need to do aspiration during the clinical injection; Long press the  For 2-5seconds ,the push rod will return 0.1ml,At the same time, observe whether there is blood return in the anesthetic bottle through the transparent booster tube. If the aspiration function is not used for a long time, it is recommended to remove the O-ring with a dental probe.

10.**Charging method**: Connect the power adapter to the charging hole at the rear of the main machine. The battery on the main machine screen shows that it is charging. At the same time, the blue indicator light on the back of the main unit lights up. It takes around 3 hours to finish charging. When the battery is full charged, the main unit screen battery shows three lights. The blue LED is off. It is recommended to charge the battery when it has one light left, which will extend the battery life.

Working system: continuous operation

No signal input and output , No anti-defibrillation function

Network power isolation: Unplug power adapter (DC 5V, 1A) for network power isolation.

Water proof grade: IPX-0

- Applied parts specified: Disposable Dental Needle connected with me equipment is chosen according doctor's needs.
- Working environment: the equipment is applied for general environment and locations .it is used at an altitude under 3000meters above with temperature of 5-40°C,relative humidity of 30-80%.
- Power adapter input: AC 100-240V , 50/60Hz Output: DC 5V, 1A.
- Input Power : 10VA
- Working power: DC 3.7V, 480mAh ,Rechargeable lithium battery
- Working voltage: DC 3.7V.
- Working current: 800mA MAX
- Stand-by time: 30 days.
- Working hours: 4hours.
- LED display time: 60 seconds
- Working instructions: LED blue light indication
- Alarm indication: LED red light indication and alarming of bi.bi
- Size: Main machine: 178mm×28mm×32mm
- Weight: Main machine: about 120g.


VII、 Operation Instructions:

1.Insert the 1.8ml card-type anesthetics (Primacaine,Scondonest.) into the booster tube according to the clinical needs, like the pic 1 .




Pic 1

2.Insert the booster tube into the fast plug of the machine, Rotate 45 degrees clockwise, like the pic 2.




Pic 2

3.Install a Disposable Dental Needle,like the pic 3.





Pic 3

4.After finish injection, turn the booster counterclockwise and pull it out, like the pic 4.



Pic 4

5.Press and hold  2-5seconds to turn on/off the power, turn the power on, the host enter standby state. double click  to reset the push rod.





6.Press and hold  2-5seconds to set the mute mode、 warning tone (beep every 0.1ml) and melody; L light flash once is mute mode, M light flash once is warning tone, H light flash once is melody; Press  once to set the speed: H light on for high speed 、 M light on for medium speed、 L light on for low speed. like the table 2.

Table 2

Speed	1.8ml Inject time
H	108s
M	150s
L	310s

7. press the button  to turn the inductive switch on/off. when the inductive switch on, the blue work indicator lights up and enters the working state; use the finger to block the inductive switch to start inject medicine, and the finger leaves the inductive switch to stop pushing the medicine .when the inductive switch off, the machine will stop inject the medicine. Press and hold the  button for 2-5 seconds in working state to do aspiration

8.Music Function: Built-in music (8M memory can be connected to the data cable to download MP3 format music or record MP3 format caring language, etc.)

Change music method:

After the main machine is connected to the computer through the Micro USB cable, you can find the internal music storage space. It contains four music files: “001 key .mp3” and “002 prompt.mp3” " 003 music.mp3" “ 004 alarm .mp3”., if you want to change the “ 003 music.mp3”, you can rename the mp3 music file name to“ 003music.mp3”, and then replace the original 003music.mp3